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Instructions for use

Intensiv ZirconCut Prosthetics: Rotating diamond-coated instruments for zirconia ceramic cutting

Prof. Mutlu Özcan, University of Zurich

Product description

- Stainless steel instruments with multiple diamond coating, placed in a stainless-steel tray for storage and protection
- Diamond grit-size: range of different sizes from coarse to fine
- Shapes: cylinder edged 90°, cylinder round, cylinder flame
- Available in FG (Friction Grip)
- Colour: black (intensity may change)
- Sterilizable and reusable

Intended use

Abrasion of zirconia ceramic and other ceramic materials in restorative Dentistry

Application details

Separation of crowns and bridges made of zirconia ceramic or other ceramic materials

Instructions for use

- Application by dentist only.
- Insert the instrument as deeply as possible into the chuck (FG) (figure a).
- Set the contra-angle at the required working speed up to a maximum value of 200'000 rpm before applying the instrument on the treated area.
- Activate abundant water spray during the entire use (minimum 50 mL/min) for cooling and rinsing away ceramic residues from the instruments.
- Ensure that there is consistent water flow, regulate in case of water shortage.
- Apply intermittently, to allow for a more thorough rinsing of the instruments and removal of ceramic residues.
- Use the instruments applying a contact force up to a maximum value of 6 N (600 g).
- Use the entire length of the working part in parallel to the surface and do not use just the tip, to allow for homogeneous wearing out of the instrument.
- It is recommended to use a dental dam for better view and safety.
- During usage periodically check that the diamond coating is uniform, intact and sufficiently present on the instrument.
- Once the application has been completed, remove the instrument from the treated area and bring it to a standstill.
- It is recommended to wear gloves and safety glasses.
- The instrument is not indicated to create endo access.

Maintenance and sterilization

- Instruments are packed non-sterile. They must be disinfected and sterilized prior to first use on the patient and disinfected, cleaned with soft brush or sonic bath and sterilized immediately after each use.
- Pre-clean by disinfection the instruments immediately after each usage.
- Disinfect the diamond-coated instruments separately from other instruments such as polymeric polishers and abrasives stones.
- Use cleaning and disinfection solutions strictly observing the concentrations and reaction times indicated by the manufacturer.
- Clean the instruments and remove possible debris after each use (with soft brush, cleaning rubber Intensiv Diakleen or ultrasonic bath), to maintain their abrasive properties.
- In case of the presence of grinding residues on the instruments it is advisable to use an ultrasonic bath, cleaning rubber Intensiv Diakleen and/or soft brush for cleaning.
- After disinfection and cleaning, inspect the instruments for residues. If necessary, repeat the disinfection/cleaning procedure.
- Cleaning and disinfection can be also carried out using professional thermosdisinfectors.

- Check for possible damage; dispose of oxidized, eccentric or deformed instruments, and instruments where the diamond-coated working part is totally or partially worn-out.
- Hot air sterilizers are not suitable for diamondcoated instruments.
- Sterilization must be carried out using the listed sterilization procedure:
 - fractionated vacuum/dynamic air removal procedure (with sufficient product drying)
 - steam sterilizer according to EN 13060/EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
 - validated according to EN ISO 17665 (valid IQ/OQ (commissioning) and product specific performance qualification (PQ))
 - maximum sterilization temperature 138 °C (280 °F; plus tolerance according to EN ISO 17665)
 - sterilization time (exposure time at the sterilization temperature) in table 1
- Storage method according to national storage indications for medical device dental instruments.
- Full cleaning, disinfection, sterilization, checking and maintenance instructions can be requested at info@intensiv.ch or visualized through the QR-code on the back page.

Risk warnings

- Always apply abundant water spray during usage and make sure that the tip of the instrument is not enclosed in the ceramic and is freely accessible by the water spray, to avoid damaging the diamond coating.
- Avoid torque, jamming or levering actions when rotating, as this increases the risk of instrument breakage.
- Never exceed the specified maximum speed, to avoid instrument breakage caused by the generation of powerful centrifugal forces.
- Partial consumption of the diamond-coated working part can cause overheating during usage of these instruments.
- Diamond coating can wear out completely during usage, and this can cause overheating. Pay attention to the state of the diamond coating throughout entire use.
- Avoid applying forces greater than the recommended values, as this could cause damage to the instrument and the treated area.
- To avoid loss of traceability of the instruments, during their entire application it is necessary to keep the outer packaging with indicated lot number.

In the unlikely case of an incident caused by the instrument please refer to the manufacturer and competent authorities.

Intensiv SA declines any liability for any use of the instrument other than that specified in these instructions for use.

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