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



Filtek = Composite; ACTIVA = Bioactive Composite; Ketac Nano = RMGI; Fuji IX = Glass Ionomer  
Source: University testing<sup>17</sup> (see back page for trademark information)

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ACTIVA = Bioactive Base/Liner; FUJI Lining = RMGI; Vitrebond = RMGI; TheraCal = Resin-Modified Calcium Silicate  
Source: Pulpdent testing<sup>18</sup> (see back page for trademark information)

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

Fig 1: Illustration shows 3-point bend test

Fig 2: Filtek = Composite; ACTIVA = BioACTIVE Composite; Ketac Nano = RMGI; Fuji IX = GI  
Source: University testing<sup>5,17</sup> (see back page for trademark information)

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



Source: Pulpdent testing<sup>20</sup> (see back page for trademark information)

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

Source: University testing<sup>1</sup> (see back page for trademark information)

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

Source: University testing<sup>16</sup> (see back page for trademark information)

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## Riferimenti

1. Fluoride ion release and recharge over time in three restoratives. University testing. Accepted for publication. Accepted for publication AADR 2014.
2. Zmener O, Pameijer CH, Hernandez S. Resistance against bacterial leakage of four luting agents used for cementation of complete cast crowns. *Am J Dent* 2014;27(1):51-55.
3. Zmener O, Pameijer CH, et al. Marginal bacterial leakage in class I cavities filled with a new resin-modified glass ionomer restorative material. 2013.
4. Flexural strength and fatigue of new Activa RMGIs. University testing. Accepted for publication AADR 2014.
5. Deflection at break of restorative materials. Unpublished university testing.
6. McCabe JF, et al. Smart Materials in Dentistry. *Aust Dent J* 2011;56 Suppl 1:3-10.
7. Cannon M, et al. Pilot study to measure fluoride ion penetration of hydrophilic sealant. AADR Annual Meeting 2010.
8. Water absorption properties of four resin-modified glass ionomer base/liner materials. (Pulpdent)
9. pH dependence on the phosphate release of Activa ionic materials. (Pulpdent)
10. Kane B, et al. Sealant adaptation and penetration into occlusal fissures. *Am J Dent* 2009;22(2):89-91.
11. Rusin RP, et al. Ion release from a new protective coating. AADR Annual Meeting 2011.
12. Sharma S, Kugel G, et al. Comparison of antimicrobial properties of sealants and amalgam. IADR Annual Meeting 2008.
13. Naorunroj S, et al. Antibacterial surface properties of fluoride-containing resin-based sealants *J Dent* 2010.
14. Prabhakar AR, et al. Comparative evaluation of the length of resin tags, viscosity and microleakage of pit and fissure sealants – an in vitro scanning electron microscope study. *Contemp Clin Dent* 2011;2(4):324-30.
15. Pameijer CH. Microleakage of four experimental resin modified glass ionomer restorative materials. April 2011.
16. Microleakage of dental bulk fill, conventional and self-adhesive composites. University testing. Accepted for publication AADR 2014
17. Comparison of Mechanical Properties of Dental Restorative Material. University testing. Accepted for publication. AADR 2014.
18. Mechanical properties of four photo-polymerizable resin-modified base/liner materials. (Pulpdent)
19. Singla R, et al. Comparative evaluation of traditional and self-priming hydrophilic resin. *J Conserv Dent* 2012;15(3):233-6.
20. Water absorption and solubility of restorative materials. (Pulpdent)
- DeRouen TA, et al. Neurobehavioral effects of dental amalgam in children: a randomized clinical trial. *JAMA* 2006;295(15):1784-1792.
- Flaim GM, et al. Remineralization of dentin lesions from a whisker-reinforced, resin-based composite. AADR 2009.
- Giorgievska E, et al. Marginal adaptation and performance of bioactive dental restorative materials in deciduous and young permanent teeth. *J Appl Oral Sci* 2008;16(1):1-6.
- Goldstep F. Proactive intervention dentistry: a model for oral care through life. *Compend Contin Educ Dent* 2012;33(6):398-402.
- Khoroushi M, Keshani F. A review of glass ionomers: from conventional glass-ionomer to bioactive glass-ionomer. *Dent Res J* 2013;10(4):411-420.
- Murray PE, et al. Analysis of pulpal reactions to restorative procedures, materials, pulp capping, and future therapies. *Crit Rev Oral Biol Med* 2002;13:509.
- Niu L, Pashley DH, Tay FR, et al. Biomimetic remineralization of dentin. *Dent Mat* 2014;30(1):77-96.
- Nordbo H, et al. Saucer-shaped cavity preparations for posterior approximal resin composite restorations: observations up to 10 years. *Quintessence Int* 1998;29(1):5-11.
- Pameijer CH, Zmerner O, Kokubu G, Grana D. Biocompatibility of four experimental formulations in subcutaneous connective tissue of rats. 2011.
- Pameijer CH, Zmerner O. Histopathological Evaluation of an RMGI ionic-cement [Pulpdent Activa], auto and light cured – A subhuman primate study. March 2011.
- Pameijer CH. Report on the retention of Embrace WetBond cement and a RMGI cement (Pulpdent). August 2012.
- Pashley DH, et al. State of the art etch-and-rinse adhesives. *Dent Mater* 2011;27(1):10.
- Peumans M, et al. Clinical effectiveness of contemporary adhesives: a systematic review of current clinical trials. *Dent Mat* 2005;21:864-881.
- Skartveit L, et al. In vivo fluoride uptake in enamel and dentin from fluoride-containing materials. *J Dent Child* 1990; 57(2):97-100.
- Spencer P, et al. Adhesive/dentin interface: the weak link in the composite restoration. *Am Biomed Eng* 2010;38(6):1989-2003.
- Spencer P, et al. Interfacial chemistry of moisture-aged class II composite restorations. *J Biomed Mater Res* 2006;77(2):234-240.
- Wang Z. Dentin remineralization induced by innovative calcium phosphate/silicate materials. AADR 2013.
- Watson TF, et al. Present and future glass ionomers and calcium-silicate cements as bioactive materials in dentistry: biophotonics-based interfacial analyses in health and disease. *Dent Mat* 2014;30(1):50-61.
- www.nidrc.nih.gov
- Yang B, et al. Remineralization of natural dentin caries with one experimental composite resin. AADR 2009.

## Letteratura Collegata

- Armstrong SR, et al. Resin-dentin interfacial ultrastructure and microtensile dentin bond strength after five-year water storage. *Oper Dent* 2004;29(6):705-12.
- Bertassoni LE, et al. Functional remineralization of dentin: induced mineral re-growth for biomechanical recovery. AADR 2009.
- Cannon ML, Comisi JC. Bioactive and therapeutic preventive approach to dental pit and fissure sealants. *Compendium* 2013;34(8):642-645.
- Comisi CC. Bioactive materials support proactive dental care. *Cosmetic Dent* 2012;1:7-13
- Delaviz Y, Finer Y, Santerre JP. Biodegradation of resin composites and adhesives by oral bacteria and saliva: a rationale for new material designs that consider the clinical environment and treatment challenges. *Dent Mat* 2014;30(1):16-32.

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